

Innovative  
and universal



## New Shaver Tools

for the **PowerDrive ART1** Motor System


## New Shaver Tools

for the PowerDrive ART1 Motor System

### New and proven




- **Sterile disposable tools**  
now also for hip arthroscopy
- **Universal resectors**  
for removing soft tissue and bone
- **Aggressive resectors**  
for arthroscopy in small joints
- **Optimized cutting designs**  
for smooth and serrated resectors
- **Disposable and reusable**  
shaver blades available
- **Automatic tool recognition**
- Preselected settings for **optimum application parameters**
- **Memory function:**  
Storage of the number of plug-in cycles and total run time



### Innovative and universal





Alongside new rotary blades and burrs Richard Wolf now supplies new **universal resectors**, which can be used for resection of soft tissue and removal of bony structures.

This system avoids an additional tool in interventions that require both soft tissue and bony structures.

Universal resectors and acromionizer for hip arthroscopy		for handles	∅ [mm]	WL [mm]	Pack [units]	Type
	Resector universal, sterile disposable	M5	4.5	160	5	499732344
	Resector universal, curved, sterile disposable	M5	4.5	160	5	499732274
	Acromionizer, sterile disposable, distal protector 2 mm reset	M5	5.5	160	5	499732565

Resectors for Arthroscopy in small joints		for handles	∅ [mm]	WL [mm]	Pack [units]	Type
	Resector, serrated reusable	M5	2.5	80	1	899750132
	Resector, serrated reusable	S1	2.5	80	1	899710132
	Resector, aggressive reusable	M5	2.5	80	1	899750332
	Resector, aggressive, sterile disposable	M5	2.5	80	5	499750332
	Resector, aggressive reusable	S1	2.5	80	1	899710332
	Resector, aggressive, sterile disposable	S1	2.5	80	5	499710332

Resectors for knee and / or shoulder arthroscopy		for handles	∅ [mm]	WL [mm]	Pack [units]	Type
	Resector, smooth reusable	M5	3.5	130	1	899750003
	Resector, smooth reusable	M5	4.5	130	1	899750004
	Resector, smooth reusable	M5	5.5	130	1	899750005
	Resector, aggressive reusable	M5	5.5	130	1	899750335



# New shaver tools

for the **PowerDrive ART1** motor system

Acromionizers for shoulder arthroscopy		for handles	∅ [mm]	WL [mm]	Pack [units]	Type
	Acromionizer reusable	M5	4.5	130	1	<b>899750554</b>
	Acromionizer, sterile disposable	M5	4.5	130	5	<b>499750554</b>
	Acromionizer reusable	M5	6.5	130	1	<b>899750556</b>
	Acromionizer, sterile disposable	M5	6.5	130	5	<b>499750556</b>

