

spirit of excellence

Perfection in Detail



Spreadable Video Mediastinoscope

Dismantlable for optimum reprocessing



spirit of excellence

Spreadable Video Mediastinoscope

Dismantlable for optimum reprocessing

Richard Wolf offers the world's first dismantlable solution in the market with the new spreadable Video Mediastinoscope. A simple twist of the wrist enables this instrument to be dismantled into its constituent elements for optimum reprocessing. Complemented by the elegant haptic design and the integrated endoscope, this mediastinoscope offers users and hygiene staff optimal prerequisites.

Extremely flexible

The spatula blades are height adjustable (parallel) and spreadable at the distal end guaranteeing optimal spatial conditions in the operating area

- Permanent optical visualization
 Continuous optical visualization
- For "expanded mediastinoscopy" The space gained by the adjustable spatula blades permits improved access to lymph nodes located laterally
- Clear view

Suction and irrigation channel for cleaning the endoscope window during the intervention

 Reduced weight and slimline design Connectable

Can connect to existing camera systems in the market

- Very easy to handle Through the use of optimized adjustment mechanism with precision ball bearings
- Can be completely dismantled Ready for reprocessing in a few simple movements of the hand
- High-quality materials All materials used are compatible with the standard commercially available reprocessing procedures, such as alkaline reprocessing, steam sterilization, sterilization with hydrogen peroxide, etc.
- In conformity with the latest reprocessing guidelines

Spreadable Video Mediastinoscope

comprising: Video Mediastinoscope with spatula spreadable distally and in parallel, WL 160 mm (8783.412), suction irrigation tube (8783.721) and tool key (9500119) 8783.4121

Additionally this purpose: Smoke evacuation tube8783.701

for fixation to the operating table
(standard rail):
U-shaped retaining arm84400.1700
Terminal socket insulated8840.9722
Fork mount

Tel.: +49 7043 35-0 Fax: +49 7043 35-4300 info@richard-wolf.com

Co-CEOs: Jürgen Pfab Jürgen Steinbeck Trade Register: Mannheim HRB 510031 VAT no.: DE144521586 Tax no.: 48020/00171

