


A clean item




## Reprocessing Baskets



to safely hold a telescope during transport,  
storage, machine reprocessing and sterilization

# Reprocessing Baskets

to safely hold a telescope during transport, storage, machine reprocessing and sterilization



- **3 different lengths**  
for virtually all endoscopic telescopes
- **Transport**  
The resilient silicone support strips protect the telescope perfectly against vibration and impacts to the basket
- **Minimal support contact area**  
Point support of telescope provides the conditions for the best possible cleaning



## Reprocessing baskets

for machine reprocessing

and sterilization,

Width: 59 mm, Height: 54 mm,

Length:

287 mm, 255 g ..... 38020.111

471 mm, 355 g ..... 38021.111

610 mm, 410 g ..... 38022.111

## Ball fixers

for reprocessing baskets

UE = 20 pcs ..... 38015.001

Technical specifications subject to change without notice.

J.610.VI.18.en.2

